

# Horizontal cylindrical pressure steam sterilizer

## Model RAU760



**East Medical Technology Co.,Ltd**

**Add:190# longping Road,Huate Industrial Zone,longgang District,Shenzhen ,GuangDong China  
518116**

**Email: info@east-medical.com**

**Tel/Fax: +86 55 33275355**

**www.east-medical.com**

## Features:

- ✧ Sterilizing course automatic sterilization controlled, easy to operate.
- ✧ With drying function, suitable dressing drying.
- ✧ With over-temperature, over-pressure auto-protect device.
- ✧ The door opening mechanism cannot be operated until the pressure in the chamber is reduced to 0.027MPa. It can't be start on if the door doesn't close well.
- ✧ The safety valve will be open automatically when the inner pressure over 0.24MPa, and the steam be exhausting to the water tank.
- ✧ The power will be cut-off automatically, cut-off the water and alarm when the machine lacking water.
- ✧ The chamber of the sterilizer is made of stainless steel.

Technical data	RAU-760	RAU-760	RAU-760	RAU-760	RAU-760	RAU-760
Sterilizing chamber	90L	150L	200L	280L	400L	500L
volume	φ440×700mm	φ440×1000mm	φ515×1000mm	φ600×1000mm	φ700×1100mm	φ700×1300mm
working pressure	0.22MPa					
working temperature	134°C					
Adjustment of	105°C - 134°C					
Time of sterilization	0-60 min					
Time of drying	0-60 min					
Heat average	≤±1°C					
Power	9KW/380V 50Hz	9KW/380V 50Hz	9KW/380V 50Hz	12KW/380V 50Hz	18KW/380V 50Hz	18KW/380V 50Hz
Dimension	1100×600×1300	1400×600×1300	1400×670×1650	1400×770×1780	1430×880×1830	1800×900×1820
Transport dimension	1250×750×1850	1550×750×1850	1560×820×1850	1680×920×2100	1800×1050×2100	1850×1050×2100
G.W/N.W	270/200 KG	320/240 KG	350/260 KG	465/365 KG	530/520 KG	580/470 KG

