



Auto CPAP DS-6

✓ Comfort

EPR(Expiration Pressure Relief,0-3 level)ensure comfort degree Automatic Synchronization Technology ensure the ventilator and patient to achieve perfect synchronization, which make patients nice experience.

✓ Quiet : Noise≤28 dB

High efficiency brushless DC motor, 10MS of rapid response Patent for cancel muffling technology.

✓ Easy Operation

Special water tank to make operation very easy and humanized.

✓ Suitable for both OSA and CSA

Auto-cpap mode can work for both OSA and CSA.

✓ Professional doctor software

Intelligent doctor software greatly reduce doctors's workload.



Specification

Description	DS-6
Work Mode	CPAP, Auto
Color	Silver
Sound	<28 dB
Screen	2.8 inch TFT, 320 x 240
Pressure Range	4-20 cmH ₂ O, auto adjustment
Data Storage	SD Card
Power Supply	Input 100-240V, 50-60Hz, 1-2A; Output + 24V, 2.5A
DC Voltage	24VDC, 2.5A
Over-Temperature Protection	Stop automatically, alarm with 4 times
Lock Function	Yes
Screen Auto Sleep	Auto sleep in 15 minutes, press any button to highlight screen
Ramp Item	0-60 min adjustable 1min/Step
Auto ON/Off	Yes
EPR	Yes, 0-3 level
Auto-Altitude	Yes
Leak Compensation	Yes
Awake Detection	Yes
Breath Event Detection	Yes
Auto Pressure Titration	Yes
Water Capacity	220ML max
Humidifier Level	Integrated Heated, 0-5 level
Alert	Power off, mask off
Dimensions(cm)	28L x 14W x 9.5H(machine); 38.5L x 36.5W x 16H(carton)
Weight(with humidifier)	1.6Kg(net weight); 3.9Kg(gross weight)
Work For	CSA, OSA