



GE Biomed

Medical System

BM 50

Biological Microscope



Product Specification:

10x plan eyepiece	- WF10x/1.8 f=25
4x achromatic objective	- NA=0, 10 WD=35.8
10x achromatic objective	- NA=0, 25 WD=6.64
40x achromatic objective	- {S}NA=0.65, WD=0.62
100x achromatic objective	- {S}NA=1.25, WD=0.35
Binocular	
Body	- Coaxial Coarse & Fine Focus Mechanism
Stage	- Size:130x140 move, Range: 78x30
Condenser	
LED parts	- 3W LED



Global Medical Engineering (BD) Ltd.

Corporate Office:17/2, Topkhana Road (2nd Floor), Dhaka - 1000, Bangladesh.

Phone: +88 09678020555 **Fax:** +88 02 9576881 **Cell:** +8801404003500

E-mail: info@gmebd.com ; gmebd@hotmail.com **Website:** www.gmebd.com

Email: info@ge-biomed.co.kr

Phone: +82(0)1067606469

Address: Incheon city, Namdong gu Gansok 3 dong,
Gansok L.H.1 danji, 105 dong 1301, South Korea