



GE Biomed

Medical System

BM 80

Biological Microscope



Product Specification:

Optical System:	Universal infinity anti-fungus optical system
Viewing Head:	30° inclined, 360° rotatable, Interpapillary distance 48 - 75 mm
Nosepiece:	Quadruple ball bearing nosepiece with positive click stops
Eyepiece:	Wide field EWF 10X /18 mm
Objectives:	BOECO Infinite corrected Plan Achromatic Objectives 4x/0.10, 10x/0.25, 40x/0.65 (S), 100x/1.25 (S, Oil)
Stage:	Large sized double layer mechanical specimen stage. 132 x 142 mm, Moving range 75 x 40 mm,
Focusing:	Coaxial coarse and calibrated fine focus control, fine division 0,001 mm, Coarse stroke 37,7mm/rotation fine stroke 0,2mm/rotation, mowing Range 20 mm safety autofocus stop unit.
Condenser:	Abbe condenser, N.A. 1.25 bright field
Illumination:	1W LED lamp, light intensity adjustable. Life time minimum 20000 hours

Email: info@ge-biomed.co.kr

Phone: +82(0)1067606469

Address: Incheon city, Namdong gu Gansok 3 dong,
Gansok L.H.1 danji, 105 dong 1301, South Korea



Global Medical Engineering (BD) Ltd.

Corporate Office: 17/2, Topkhana Road (2nd Floor), Dhaka - 1000, Bangladesh.

Phone: +88 09678020555 **Fax:** +88 02 9576881 **Cell:** +8801404003500

E-mail: info@gmebd.com ; gmebd@hotmail.com **Website:** www.gmebd.com