

▶ Test Menu

▶ Thyroid

- TSH
- T3
- T4
- FT3
- FT4
- Anti-Tg*
- Anti-TPO*

▶ Allergy

- IgE

▶ Infectious Disease

- TB-IGRA
- HBs-Ag
- Anti-HCV
- HIV Combo
- TP

▶ TORCH

- TOXO IgG*
- TOXO IgM*
- RUB IgG*
- RUB IgM*

▶ Fertility

- β -HCG
- FSH
- LH
- PROL
- PROG
- TESTO
- E2
- AMH
- E3*
- INH-B*

▶ Brain Injury

- S100
- A β 1-42*
- AD7C-NTP*

▶ Bone Metabolism

- 25-OH-VD
- PTH
- Calcitonin*
- N-MID Osteocalcin*

▶ Tumor Markers

- CEA
- AFP
- tPSA
- fPSA
- CA125
- PG I
- PG II
- Gastrin-17
- CA19-9
- CA15-3*

▶ Anemia

- VB12
- Folate
- Ferritin

▶ Inflammation

- PCT
- CRP
- SAA
- IL-6
- HBP

▶ Cardiac Markers

- cTnI
- Hs-cTnI
- cTnT
- Myo
- CK-MB
- NT-ProBNP
- D-Dimer
- H-FABP
- BNP
- Lp-PLA2*
- ST2*
- MPO*
- Hcy*
- GDF-15*

▶ Diabetes

- Insulin
- C-Peptide

* Upcoming

Global Medical Engineering (BD) Ltd.

Corporate Office: 17/2, Topkhana Road (2nd Floor), Dhaka - 1000, Bangladesh.

Phone: +88 09678020555 Fax: +88 02 9576881 Cell: +8801404003500

E-mail: info@gmebd.com ; gmebd@hotmail.com Website: www.gmebd.com

All rights reserved. Features and specifications are subject to change without prior notice.



AIKANG



CLIA-mate 1200

Chemiluminescence Immunoassay Analyzer

Aikang MedTech Co., Ltd



**2 Chambers, 12 Channels
Free Combination of Test Items**

Output

- ⊙ Built-in printer
- ⊙ Available for the external printer
- ⊙ Available for LIS

Super Touchscreen

- ⊙ 15.6-inch color touchscreen
- ⊙ Show all 12 testing in one screen
- ⊙ User-friendly interface design
- ⊙ Simple and concise software

Integrated Reagent Chamber

- ⊙ Two independent chambers, 12 testing channels in total
- ⊙ Highly integrated reaction module, no liquid path
- ⊙ Automated reagent barcode scanning
- ⊙ Disposable tip for each assay
- ⊙ Allow more different items running simultaneously
- ⊙ 15 minutes for the first result


Compact

- ⊙ Dimensions: 487*573*631 mm
- ⊙ Weight: 60 kg


Single Test Reagent

- ⊙ Specifications: 30 test/kit, 60 test/kit
- ⊙ Ready to use, no worry of onboard stability
- ⊙ Easy reagent registration by barcode scanning
- ⊙ Calibration validity 100 days, calibrator & QC included in the kit package
- ⊙ Comprehensive test menu

All-in-one Reagent Cartridge



Easy Operation



Step 1: Sample adding



Step 2: Reagent cartridges loading



Step 3: Button start