


## L6 LED Surgical Light



### Physical Specifications

L6 LED Surgical Light	
Dimension	Width: 358 cm (level cross arm, lamp panel, display, and camera are all in a vertical direction) Height: 171 cm (level cross arm, lamp panel, display, and camera are all in a vertical direction) Length: 473 cm (level cross arm, lamp panel, display, and camera are all in a vertical direction)
Weight	85kg
Control Panel	
Dimension	140 mm×127 mm×36.8 mm
Touch Screen	
Dimension	4.3 inches (display area: 95 mm × 54 mm)
Resolution	480×272 pixels

### Shenzhen Comen Medical Instruments Co.,Ltd

HeadquarterAdd:Floor7,Block5,  
4th Industrial Park of Nanyou,Nanshan District,Shenzhen City 518052,China.  
E-mail : [service@szcomen.com](mailto:service@szcomen.com)

Brightness	Non-adjustable		
Range of Rotation	60° ( $\pm 2^\circ$ )		
<b>Rotation Angle and Startup Force of the Equipment</b>			
	<b>Items</b>	<b>Rotation Angle</b>	<b>Startup Force</b>
<b>Mother light</b>	Cross arm around the column	Greater than 360°	Less than 35 N
	Spring arm around the cross arm	340° ( $\pm 5^\circ$ )	Less than 35 N
	Spring arm moving up and down	45° (upper), 50° ( $\pm 5^\circ$ ) (lower)	Less than 35 N
	Upper elbow around the spring arm (X direction)	340° ( $\pm 5^\circ$ )	Less than 15 N
	Lower elbow around the upper elbow (Y direction)	340° ( $\pm 5^\circ$ )	Less than 15 N
	Lamp panel around the lower elbow (Z direction)	180° ( $\pm 5^\circ$ )	Less than 15 N
<b>Son light</b>	Cross arm around the column	Greater than 360°	Less than 35 N
	Spring arm around the cross arm	340° ( $\pm 5^\circ$ )	Less than 35 N
	Spring arm moving up and down	45° (upper), 50° ( $\pm 5^\circ$ ) (lower)	Less than 35 N
	Upper elbow around the spring arm (X direction)	340° ( $\pm 5^\circ$ )	Less than 15N

**Shenzhen Comen Medical Instruments Co.,Ltd**  
HeadquarterAdd:Floor7,Block5,  
4th Industrial Park of Nanyou,Nanshan District,Shenzhen City 518052,China.  
E-mail : [service@szcomen.com](mailto:service@szcomen.com)