



ZEUS-400P Surgical unit

Most Safe Realiability



Description

Feature

ZEUS-400P is the up-graded version of ZEUS-400P which is a stat of the art technology electrosurgical unit for universal purpose.

-5 kinds of cut, bipolar and monopolar coagulation.

-Best safty with Feedback system

REM(Return Electrode Monitoring)monitors the size of contacting area between a patient and the pad, if the size is inappropriate, it automatically blocks the high-frequecy current to minimize the danger of burning incidents.

-Stable output power

Unit can stably and stagithly output power by using micro processor

When the unit turn it off and turn it back on, last set value will be displayed by flash memory.

-Available bipolar cut of 100W

-compatible with Argon-Z (Argon Plazma Coagulation)

ZEUS-400P Surgical unit

Model		ZEUS-400P
CUT	Pure cut	400W / 300Ω
	Blend 1	250W / 300Ω
	Blend 2	200W / 300Ω
	Blend 3	150W / 300Ω
	Contact	120W / 300Ω
Coagulation	Spray	100W / 300Ω
		100W / 100Ω
Bipolar Coagulation		100W / 200Ω
Bipolar Cut		
Main Frequency	400, 500 KHz	
Power Source	AC120V or AC 230V, 50 / 60 Hz	
Protection Class	Class 1, Type CF	
Dimension(WxDxH)	330x455x150 mm	



3 Pin plug cord/cable & Removable stainless Steel Ultra-clean Blade Electrode Safe Use, flexible, and economical



Bipolar Forceps

The standard coagulation forceps for the use of safe bipolar technology have a voltage-resistant safety coating



Patient-return-plate-cable-L-3M-plug

ZEUS connecting cables for electrosurgery are designed for the high demands of everyday surgical procedures. They provide a high-level of electric strength, safety stripes, and has a high reusability.



Bipolar Cable

Safe, flexible, and economical



SS Patient Plate

The reusable neutral electrodes made of conductive SS close the monopolar circuit on the passive side between the electrosurgical unit and the patient.



"Foot switch"

Foot switches are used to activate electrosurgical instruments and are connected to ZEUS

Manufacturer:

Zerone Co. Ltd

13 LS-Ro, Gunpo-si, Gyeonggi-do,
City Gyeonggi-do, Korea.
Phone: 82 - 31 - 6895621,
Fax: 82 - 31 - 4535444
Web: www.01zeus.com



Local Distributor:

Global Medical Engineering (BD) Ltd.
17/2 (2nd Floor), Topkhana Road, Dhaka-1000, Bangladesh
Phone: +88 02 9564225, Fax: +88 02 9576881
E-mail: info@gmebd.com, gmebd@gotmail.com
Web: www.gmebd.com