



NEW

DP-20

Digital Ultrasonic Diagnostic Imaging System

When Expectations Become Reality

mindray
healthcare within reach

DP-20

Digital Ultrasonic Diagnostic Imaging System

DP-20, a smart portable B/W ultrasound system from Mindray, brings more advances to general practices to improve patient's benefits. With deep understanding of needs from customers, we hold a belief that this new generation of portable ultrasound system will be your best choice when you want more while cost less.



Brilliant Ergonomics

- High resolution 12.1" LED with tilting functionality
- User-friendly keyboard and controls

Excellent Image Quality

- Tissue Harmonic Imaging (THI) enhancing contrast resolution
- Quick image optimization by IP (Image Processing)
- 8-segment TGC allowing delicate image adjustment


Powerful Workflow

- iStation, One-Key-Save, user-defined keys bringing great efficiency in your daily diagnosis

Clinical Versatility

- 8 probes with abundant measurement packages covering traditional ultrasound applications and emerging fields such as urology, MSK and anesthesia


PROMOTE




Median Nerve



Thyroid



Hepatic Hemangioma



Carotid

REDUCE

Weight

- Only 5.5kg. Take it to your patients.

Cost

- Diagnosis needs and budget conflict? Have a try on Mindray DP-20.

DISTRIBUTOR:



MINDRAY is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Specifications subject to changes without prior notice.
© 2010-2011 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.
P/N:ENG-DP-20-210285-20110901

mindray

Mindray is listed on the NYSE under the symbol "MR"
Mindray Building, Keji 12th Road South, High-tech Industrial Park,
Nanshan, Shenzhen 518057, P.R. China
Tel: +86 755 8188 8998 Fax: +86 755 26582680
E-mail: intl-market@mindray.com Website: www.mindray.com