



**DP-5** **NEW**

Digital Ultrasonic Diagnostic Imaging System

**Ultimate Master**



**mindray**

healthcare within reach

# DP-5

## Digital Ultrasonic Diagnostic Imaging System

### Guaranteed Diagnostic Confidence

#### iClear™ Speckle Suppression

Improve contrast resolution without compromising spatial resolution

#### iBeam™ Spatial Compounding Imaging


Significantly improve resolution and meanwhile get more information of tissue

#### Phase Shift THI

Increase noise rejection capability,

Better contrast resolution,

Better tissue subtleties visualization



## DP-5, Ultimate Master

Since 2002 when Mindray launched Chinese first B/W ultrasound system DP-9900, we have never stopped to improve our products in accordance with our customers' needs. Now we offer you the new DP-5, the Ultimate Master, which deeply reflects our sustained efforts and remarkable achievements in developing B/W ultrasound system.

Combined with the latest imaging technologies, wonderful ergonomic design and intelligent workflow, DP-5 will bring you brand new experience and become a powerful and reliable assistant in your daily diagnosis.



## Comfortable Operation Experience

High resolution, wide angle 15" LCD

Compact and streamlined shape

User friendly keyboard and controls

Light weight probes of wide range of broadband



## Great Productivity

### iTouch™ Auto Optimization

Optimizes the uniformity and dynamic range with just one button

### iZoom™

Instant Full Screen Viewing

### Quick Save

Image transfer to local or directly to USB with one push

### iStation™

integrates patient data storage, archiving, review and retrieval efficiently

# DP-5

## Digital Ultrasonic Diagnostic Imaging System

### Mindray Medical Netherlands B.V.

Drs.W.van Royenstraat 8, P.O. Box 26 3870 CA Hoevelaken,  
The Netherlands  
Voice: (31-33) 254-4911  
Fax: (31-33) 253-4280

### Mindray Medical Italy S.r.l.

Via Girardo Patecchio, 4 20141 Milano Italy  
Voice: (39-02) 57-402-444  
Fax: (39-02) 5523-2018

### Mindray Medical España S.L.

Avenida Manteras 38, Edificio B, Bajo B008. 28050-Madrid Spain  
Voice: (34-91) 392 3754  
Fax: (34-91) 392 3756

### Mindray Medical Germany GmbH

Zwischen den Bächen 4, 64625 Bensheim, Germany  
Voice: (49) 6251-17524-0  
Fax: (49) 6251-17524-20

### Mindray (UK) Limited

3 Percy Road, St Johns Park, Huntingdon, Cambs, PE29 6SZ  
United Kingdom  
Voice: (44) 1480 416 840  
Fax: (44) 1480 436 588

### Mindray Medical France SARL

Europarc Creteil 1 allée des cerisiers 94000 CRETEIL, France  
Voice: (33-1) 4513-9145  
Fax: (33-1) 4513-9151

### Mindray Medical Rus Co.Ltd

2 Zvenigorodskaya, 13, bld 41, 123022 Moscow, Russia  
Voice: (7-499) 553-6036  
Fax: (7-499) 553-6039  
Hotline: 8-800-333-5323

### Mindray Medical Colombia S.A.S

Room No.702,Torre UNIKA, Carrera 9 No.77-67, Bogota, Colombia  
Voice: (57-1) 3130892 / 3210916

### Mindray do Brasil – Comércio e Distribuição de Equipamentos Médicos Ltda.

Rua Tavares Bastos, 329, Perdizes, São Paulo - SP, CEP 05012-020  
Telephones: (55 11) 3124-8026

### Mindray Medical Mexico S. de R.L. de C.V.

Félix Parra # 175 Colonia. San José Insurgentes Sur  
Delegación Benito Juárez 03900 México, D. F.  
Voice: (52-55) 5661-9450 / (52-55) 5662-6620  
Fax: (52) 55 56 62 65 97

### PT.Mindray Medical Indonesia

Menara Dea 1,Suite 301,Kav.E4 No.3/1, Mega Kuningan,  
Jakarta Selatan 12950  
Phone: +6221 5762650  
Fax:+ 6221 5762655

### Mindray Medical Egypt

Ground Floor, 11 El Imam Metwally El Sharawi Street, Sheraton Heliopolis - El Nozha, Cairo, Egypt  
Tel : +2 02 22691047  
Fax : +2 02 22691046

### Mindray Medical India Pvt. Ltd. (Delhi Office)

401/402, NDM-1, Netaji Subhash Place, Wazirpur District Center,  
New Delhi - 110034  
Voice: +91 11 4923 0000  
Fax: +91 11 4923 0030  
Toll Free - 1800 10 20 300

### Mindray Medical India Pvt. Ltd. (Mumbai office)

B-404, City Point, Andheri Kurla Road, Andheri East, Mumbai,  
400 059, India  
Voice: (91-22) 4020-0000  
Fax: (91-22) 4020-0011

### Mindray Medikal Teknoloji Istanbul Ltd. Sti.

Yenibosna, Köyaltı Mevkii, 29 Ekim Cad. Kuyumcukent 2 Nolu Plaza Kat:2 Daire:1, Yenibosna, Istanbul, Turkey  
Voice: (90-212) 482 08 77



● Subsidiary ● R&D center



MINDRAY is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd.  
Specifications subject to changes without prior notice.  
© 2011 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.  
P/N: ENG-DP5-420285X4-20111130

# mindray

Mindray is listed on the NYSE under the symbol "MR"  
Mindray Building, Keji 12th Road South, High-tech Industrial Park,  
Nanshan, Shenzhen 518057, P.R. China  
Tel: +86 755 8188 8998 Fax: +86 755 26582680  
E-mail: intl-market@mindray.com Website: www.mindray.com