



# family cares



Become the first choice  
in respiratory care

---

# RESmart GII

## RESmart GII CPAP System



### E-20C-H-0

RESmart GII CPAP System is a basic fixed-pressure device with independent humidifier, high heating efficiency and easy data access by BMCares iCode.

## RESmart GII Auto CPAP System



### E-20A-H-0/E-20AJ-H-0

RESmart GII Auto CPAP System is a premium auto-adjusting pressure device with independent humidifier, high heating efficiency and easy data access by BMCares iCode.

## RESmart GII BPAP System



**RESmart GII S** BPAP System is ideal for sleep apnea patients who find it difficult to adjust to a continuous fixed-pressure machine.

**RESmart GII AutoS** BPAP System is a premium auto-adjusting pressure device providing enhanced pressure by advanced technologies.

**RESmart GII ST** BPAP System provides non-invasive ventilation with complete ventilation modes and volume support function to target tidal volume.

Model	Working mode	Pressure range (cmH <sub>2</sub> O)	Screen	Device type
E-20C-H-0	CPAP	4 - 20	2.4-inch	CPAP
E-20A-H-0	CPAP, AutoCPAP	4 - 20	3.5-inch	Auto CPAP
E-20AJ-H-0		4 - 20	2.4-inch	
T-20S	CPAP, S	4 - 20	2.4-inch	Bi-level S
T-25S		4 - 25	3.5-inch	
T-20A	CPAP, S, Auto S	4 - 20	2.4-inch	Auto Bi-level
T-25A		4 - 25	3.5-inch	
T-20T	CPAP, S, T, S/T	4 - 20	2.4-inch	Bi-level S/T
T-25T		4 - 25	3.5-inch	
T-30T		4 - 30	3.5-inch	
U-20T		4 - 20	2.4-inch	
Y-25T	CPAP, S, T, S/T	4 - 25	3.5-inch	Bi-level S/T with target tidal volume
Y-30T		4 - 30		
Y-30AT	CPAP, AutoCPAP, S, Auto S, T, S/T	4 - 30	3.5-inch	Full mode lab titration device



- CPAP System
- Auto CPAP System
- BPAP System
- Lab Titration Device

- ✓ Obstructive Sleep Apnea (OSA)
- ✓ Central Sleep Apnea (CSA)
- ✓ Mixed Sleep Apnea
- ✓ Periodic Breathing
- ✓ Respiratory Insufficiency

- Integrated humidifier:  
high heating efficiency
- Push-Push design:  
easy way to disassemble/ install water tank
- Built-in cellular module:  
makes connectivity standard
- Preheating system
- Integrated heating tubing
- Noise level: 28 dB
- Smart A/B/C technologies

## ASV ventilator

Designed for patients with Cheyne-Stokes respiration / periodic breathing, learning the patient's breathing state quickly, maintaining normal minute ventilation and stabilizing breathing.

Model	Pressure range (cmH <sub>2</sub> O)	Working mode					Targeted tidal volume
		CPAP	Auto CPAP	Auto S	S	T, S/T	
G3 C20	4-20	✓					
G3 A20	4-20	✓	✓				
G3 B20A	4-20	✓		✓	✓		
G3 B25S	4-25	✓			✓		
G3 B25A	4-25	✓		✓	✓		
G3 B25VT	4-25	✓			✓	✓	✓
G3 B30VT	4-30	✓			✓	✓	✓
G3 B30SV	4-30	✓			✓	✓	✓
G3 LAB	4-30	✓	✓	✓	✓	✓	✓

## RESmart G2S 2018



### RESmart G2S CPAP System G2S C20

RESmart G2S CPAP System is a basic fixed-pressure device with integrated humidifier and easy data access by BMCares iCode.

### RESmart G2S Auto CPAP System G2S A20

RESmart G2S Auto CPAP System is an auto-adjusting pressure device with integrated humidifier, high heating efficiency and easy data access by BMCares iCode.

## RESmart G2S BPAP System



### RESmart G2S S BPAP System

- G2S B20S
- G2S B25S

### RESmart G2S AutoS BPAP System

- G2S B20A
- G2S B25A

### RESmart G2S S/T BPAP System

- G2S B20T
- G2S B25T

### RESmart G2S Lab Titration Full Mode Device

- G2S LAB

## RESmart



RESmart CPAP System



RESmart Auto CPAP System



RESmart BPAP System

## BMC product other features

- SpO<sub>2</sub> monitoring
- Ramp
- ECO Smart Heating technology
- Multi-language
- Auto On
- Auto Off
- Tubing: 22mm or 15mm
- Manual P
- Delay Off
- Power failure alert
- Mask Test
- GPRS/ WIFI Connectivity
- Rise Time
- BMCares iCode data access



## BMCare's iCode Virtual is real




Advanced Clinical Efficacy Reporting with QR Coded Data for Smartphone Uploaded  
Reports Over 1, 7, 30, 90, 182, 365 Days

GENERATIONS	iCode	iCode QR	iCode QR+
Terms			
Format	Strings	QR	QR

Automatically send reports to linked emails (up to 3)

## BMCare's Cloud

Broader and Deeper, BMCare's bridge the patients and clinicians via cloud-based platform. Respiratory health management couldn't be easier. Attach the wireless module to your device, you can stay more connected to avoid unnecessary visits.



Device record and upload real time data

Diagnosis based on real time records

## Others



● SpO<sub>2</sub> Kit



● Power Bank (D1402)



● DC/DC Converter



● GPRS / Wi-Fi Kit



● Heated Tubing



● Ventilator Disinfector

## BMC Medical Co., Ltd.

Add: Room 110 Tower A Fengyu Building, No. 115 Fucheng Road, Haidian, 100036 Beijing, PEOPLE'S REPUBLIC OF CHINA  
Tel: +86-10-51663880 Fax: +86-10-51663880 ext.810 Website: [en.bmc-medical.com](http://en.bmc-medical.com) Email: [intl@bmc-medical.com](mailto:intl@bmc-medical.com)