



hemochroma PLUS

Total Hemoglobin meter

- ▶ Fast and Simple
- ▶ Cuvette with Dual Sample Inlets
- ▶ Capillary blood + Tube blood
- ▶ US FDA CLIA Waived



hemochroma PLUS is an in vitro diagnostic device dedicated to quantitative testing of total hemoglobin concentration in the blood. It is lightweight with small dimensions and battery operated. It is easy to operate and convenient. It provides accurate and reliable inspection results compared with large analyzers

Specification

Item	Specification
Cat#	FPRR005
Dimensions	106.5 mm (L) x 151.0 mm (W) x 38.5 mm (H)
Weight	215 ± 3 g
Power supply	Battery (AA 1.5V x 4ea) / Adaptor / USB power
Data output	LCD, thermal printer (option)
Operation temperature	15 ~ 35°C (59 ~ 95°F)
Measuring time	3 sec
Detection range	5.0~25.6 g/dL
Sample	Whole blood (Capillary, tube blood (15 µL))
Memory	Test results up to 1,000
Interface	RS-232 1 Port USB
Display LCD	3.2" TFT color LCD