

01

HOLTER SYSTEM



Introduction:

Holter ECG product was initiated by Mr Norman Jefferis Holter in 1949. Holter Recorders can record at least more than 24 hours ECG waveforms, including data when patients rest, exercise, dine, work, study and sleep under different conditions. Can find out Cardiac Arrhythmia and myocardial Ischemia etc heart disease much easily, is the most important basis of clinical analysis, diagnosis establish and the curative effect judgment.

Hardware



- **Model:**DMS300-4A
- **Channel number:** 12 Lead
- **Recording time:** $\approx$ 3 Days
- **Storage capacity:** $\approx$ 1G
- **Size:**8.8CM x 5.5CM x 2.1CM
- **Power supply:** 1 x AAA battery



- **Model:**DMS300-3A
- **Channel number:** 3 Lead
- **Recording time:** $\approx$ 7 Days
- **Storage capacity:** $\approx$ 1G
- **Size:**8.8CM x 5.5CM x 2.1CM
- **Power supply:** 1 x AAA battery

Mini dynamic ECG recorder - features introduction:

- Independent pacing channel, pacing channel 10000Hz;
- Normal ECG sampling rate is 4096Hz;
- 3 leads and 12 leads are identified;
- USB2.0 fast playback;
- Support DMS ECG information management system ;
- Special event button, ECG acquisition indicator;



05

- Support telephone electrocardiogram transmission;
- Support late potential data acquisition;
- Support vector ECG data acquisition;



- **Model:**DMS300-4AL / DMS300-3AL
- **Channel number:** 12Lead
- **Recording time:** $\approx$ 5 Days
- **Storage capacity:** $\approx$ 1G
- **Power supply:** 1 x AAA battery

- DMS Holter ECG recorder (DMS300-4AL) has passed the international mainstream organization testing and certification of FDA、CFDA、CE、ISO13485.
- In conjunction with supporting software,can record and analyse 12-lead 24-hour or multi-day ECG data.
- A variety of sensors and high sampling rates ensure that the data content is rich, accurate, medical and scientific institutions can then provide high-quality services.



- **Model:** DMS300-4A / DMS300-3A
- **Channel number:** 12Lead / 3Lead
- **Recording time:** 24 hours
- **Storage capacity:** $\approx$ 1G
- **Size:** 8.8CM x 5.5CM x 2.1CM
- **Power supply:** 1 x AAA battery

DMS CUS version of the Holter ECG recorder has passed the detection and certification of the international mainstream organizations (FDA、CFDA、CE、ISO13485) .In conjunction with supporting software,can record and analyse 12-lead 24-hour or multi-day ECG data.

06

Component



3-Lead ECG Cable



12-Lead ECG Cable



Data Transfer Cable



Pouch For
DMS300-3A/300-4A



Pouch For
300-3AL/300-4AL

