

HYHJ series

Video laryngoscope

Application:

Anesthesia
ICU
Emergency Center
Pediatrics
Ambulance
Intubation training
.....



Video laryngoscope

Power Button

ON/OFF&freeze

Display

3 inch HD
screen,UP/DOWN,LEFT/RIG
HT turning

Alloy handle

Ergonomic design,easy handle

Medical stainless steel

suitable for different
kind of disinfection.reusable

USB interface

charging and data transfer

Shooting

short press for shooting,
long press for recording

2 million pixels CMOS



One unit with multiple models of blades (child,adult,difficult airway)

disposable blade



reusable blade



product description

Considering clinical application,use high-class material and advanced technology,updated the design and picture quality.

product feature

simple but nice design,clear picture,easy operation,shooting,recording,HDMI,data transfer,long-life battery

Configuration

the monitor can connect to any kind of blade we supply.5 size total,for infant,pediatrics,child,adult,difficult airway.



XIAN HAIYE MEDICAL EQUIPMENT CO.,LTD.

Add: Fengjing Industrial Park, Xian City, Shaanxi Province, China. Tel/Fax: 029-85416591

Email:sale@haiyemedical.com Web:en.haiyemedical.com