

i-chamber



The i-chamber is an auxiliary device for the ichroma (ichroma™, ichroma™ II) systems and provides an environment set to the correct temperature to the ichroma cartridge when heating is required between tests. You can improve the reliability of your test results by keeping the temperature constant, one of the environmental factors that can cause unexpected errors in your test results.

Specification

Item	Specification
Cat#	FPRR009
Dimensions	190.0 mm(L) x 120.6 mm(W) x 108.0 mm(H)
Weight	1.5 kg (± 0.3 kg)
Power supply	DC 12V, 5A
Data output	LCD (16 x 4 characters)
Operation temperature	15 ~ 35°C (59 ~ 95°F)
Slots	6
Interface	Mini USB 1 Port