

Infusion Pump



- 4" large LCD display with big font can be seen clearly 5m away; Ergonomics operation interface make it easy to adjust and monitor the parameters so as to reduce
- Mistake operation.
- Flow Velocity Mode, Drop Rate Mode and Time Mode for wide range of applications.
- Micro infusion speed 0.1ml/h up to high speed 1500ml/h satisfies special clinic scenario including
- Requirements for Neonatal.
- Software adjusting function suit for any brand of infusion apparatus.
- Automatic light up upon door open for easy operation at night.
- 12 V DC power supply, nurse call interface available.
- Light weight: 1.8Kg.

Model: ME600

Origin: China

Manufacturer: Shenzhen Comen Medical Inst,