


M201 Microplate Reader



Technical Specifications

System function

Readout range	0.000~4.000Abs
Plate type	96-well plate
Measure system	8-channel optical system
Wavelength	405, 450, 492, 630nm, two more options
Band width and accuracy	8nm ± 2nm
Light source	Monochromatic LED
Reading speed	5s with single wavelength for whole plate
Repeatability	< 1%
Linearity	R ≥ 0.995
Stability	< 0.005Abs/10 min
Shaking	3 shaking modes with variable time
Input	Mouse, Membrane keyboard, Touch screen

Display

7" TFT LCD, 800 × 480 pixels

Printer

Built-in printer, external printer

Working conditions

Power supply	AC 100~240V, 50/60Hz
Temperature	10°C~35°C
Humidity	15%~85%

Interface

RS-232, 4 USB ports, SD card

Dimension (mm)

400(L) x 310(W) x 169(H)

Weight

6 kg

Options



Keyboard



External printer



Touch screen



Labfriend software



Barcode reader



Flash disk



M201
Microplate Reader

EMP
Shenzhen Emperor Electronic Technology Co.,Ltd.

Add:2/F, Building7, Tian'an Nanyou Industrial Area, Nanshan District, Shenzhen 518054, China

Tel: +86-755-26071664 26415597

Fax: +86-755-26073886 26419886

E-mail: business@china-emperor.com

Http://www.china-emperor.com



We reserve the right to make modifications without prior notification.

DISTRIBUTOR

Global Medical Engineering (BD) Ltd.
Corporate Office:17/2, Tophkhana Road (2nd Floor), Dhaka - 1000, Bangladesh.
Phone: +88 09678020555 Fax: +88 02 9576881 Cell: +8801404003500
E-mail: info@gmebd.com ; gmebd@hotmail.com Website: www.gmebd.com




Innovation with LED light source

M201 Microplate Reader

Designed with LED light technology, M201 is green, economical and in line with the trend of environment friendly development

Features

- 7" TFT LCD screen
- Test results simultaneous displayed in screen
- LED light source with lifespan over 100,000 hours
- Bichromatic, 8-channel optical system
- Large storage capacity up to 200 test protocols and 100,000 results
- 3 shaking modes with variable time
- Windows operation system
- Self-check, malfunction alarm



M201 Microplate Reader



Applications

- Pesticide residue in food
- Melamine
- Hormone, such as T3, T4 and TSH
- Test Giant cell IgG antibody, Toxoplasma IgG antibody, Rubella IgM antibody, Fever blister I-IgM antibody and Fever blister II-IgM antibody
- Cancer marker, Typhoid fever, Paratyphoid fever
- Sexually transmitted disease (Gonorrhoea, Syphilis, HIV) and Hepatitis (Hepatitis-A, Hepatitis-B, Hepatitis-C)
- Fetal antibody and vaccinum reaction, etc

