

Operation table 3008

The operation table 3008 is supplied for surgical operation of chest, abdomen, E. E. N. T., gynaecology and obstetrics, orthopaedics department.

It is elevated by foot-oil pump, mechanical transmission are applied for head operation in other parts of the unit. The feet in lay-bar is made of stainless steel.

Our operation table 3008 also have other models:

3008A: Base and the cover both made of stainless steel

3008A-1: With radiography stand, base and the cover both made of stainless steel

3008A(01): Table height adjustable by electric control, base and cover made of stainless steel

3008B: Tabletop, base and the cover both made of stainless steel

3008B-1: With radiography stand, tabletop, base and the oil pump's cover both made of stainless steel



v id="StyleTableProd">

Item	Contents
legth & wide	2100mm x 480mm
max. & min. height	1045mm - 800mm
forward & backward	> 35
left & right	> 25°
Oil pump Vertical of Oil pump	240mm
head up & down	90°
back up> 85° back down> 10°	
elevating height of kidney bridge	> 120mm
leg down	> 90°

Larkcop International Co., Ltd.

City: shanghai

Province/State:shanghai

Country/Region:China

Tel:86-21-68543228