



RESmart GII BPAP System (T Series)

Bi-level positive airway pressure with Ti Min/Ti Max

RESmart GII BPAP Series provide non-invasive ventilation for patients with Obstructive Sleep Apnea (OSA) and respiratory insufficiency in hospital or at home. Ti Min/Ti Max ensure inspiration time and improve the treatment effect significantly.



3.5-inch colorful LCD Display

Real time pressure & flow wave display.

SpO₂ Kit

Blood oxygen saturation and pulse rate display.

Eco Smart Technology

Eco Smart heating system with innovative dual water chambers design, accurate water quantity control, with real time compensation. Eco Smart warm up, high efficiency, lower water bills.

iCode QR



Compliance data is encoded by QR code. Search “BMCares” on Google Play or App Store.

Ti Min/Ti max

Ensure inspiration time, improve the treatment effect significantly.

Multiple data retrieval

SD Card / GPRS / Wi-Fi / iCode / iCode QR / iCode QR+



	T-20S	T-20A	T-20T	T-25S	T-25A	T-25T	T-30T
Mode	CPAP, S	CPAP, S, Auto S	CPAP, S, T, S/T	CPAP, S	CPAP, S, Auto S	CPAP, S, T, S/T	CPAP, S, T, S/T
Pressure		4 - 20 cm H ₂ O			4 - 25 cm H ₂ O		4 - 30 cm H ₂ O
LCD Screen		2.4 inches			3.5 inches		
Weight		5.5 lb (2.5 kg) with humidifier			5.5 lb (2.5 kg) with humidifier		
Dimensions (L × W × H)		11.4" × 7.7" × 5.3" (290 × 196 × 134 mm) with humidifier			11.4" × 7.7" × 5.3" (290 × 196 × 134 mm) with humidifier		
SpO ₂ Kit		Optional			Optional		
Data Storage		SD Card			SD Card		
Reslex		Off / 1 / 2 / 3			Off / 1 / 2 / 3		
Ramp P		4 cm H ₂ O			4 cm H ₂ O		
Max Ramp		0 - 60 minutes			0 - 60 minutes		
Ramp Time		0 - Max Ramp			0 - Max Ramp		
Humidity Level		off / 1 - 5			off / 1 - 5		
Ti Min		0.1 s - Ti Max			0.1 s - Ti Max		
Ti Max		0.3 s - 3.0 s			0.3 s - 3.0 s		

BMC Medical Co., Ltd.