

M-6

Tube Roller



Features

Speed range of 0-70rpm

Six roller design for gentle rolling and easy to clean

Tube rollers provide a gentle and highly efficient rolling motions, ideal for mixing blood samples, viscous substances and liquid-solid suspensions. Used in a variety of applications, such as the prevention of blood coagulation, immune precipitation, etc.

Specifications		M-6
Amplitude		24mm
Motor type		DC motor
Rocking motion		Rocking and rolling
Max. load capacity		4kg
Number of rollers		6
Roller size[length]		280mm
Speed range		0-70rpm
Speed display		Scale
Timer		-
Timer display		-
Time setting range		-
Operation mode		Continuous
Voltage		100-240V,50/60Hz
Power		25W
Dimension		450×260×120mm
Weight		4.5kg
Permissible ambient temperature and humidity		5-40 °C,80%RH
Protection class		IP21



Global Medical Engineering (BD) Ltd.

Corporate Office: 17/2, Topkhana Road (2nd Floor), Dhaka - 1000, Bangladesh.
Phone: +88 09678020555 **Fax:** +88 02 9576881 **Cell:** +8801404003500
E-mail: info@gmebd.com ; gmebd@hotmail.com **Website:** www.gmebd.com

GE-BIOMED

+82(0)1067606469

www.ge-biomed.co.kr

Incheon city, Namdong gu Gansok 3 dong,
 Gansok L.H.1 danji, 105 dong 1301, South Korea