

Intuitive operation interface: Configured with 3" LCD and large font for greater accessibility and visibility.

Sleek ergonomic design:
Only 1.6kg with handle



High adaptability design:
Multifunctional connector for
RS232, AC, nurse call.

• Longer Battery Life: No less than 12hrs*

*Under specific conditions.



M260 Syringe Pump



Shenzhen Comen Medical Instruments Co., Ltd.

Add: Floor 10, Floor 11 and Section C of Floor 12 of Building 1A & Floor 1 to Floor 5 of Building 2, FIYTA Timepiece Building, Nanhuan Avenue, Matian Sub-district, Guangming District, Shenzhen, Guangdong, 518106, P.R. China
Tel:+86-755-2640 8879 Fax:+86-755-2643 1232 Website: en.comen.com E-mail: info@szcomen.com

P/N: EN-M260-8P-20220621-V1.0

Our Latest Frugal Innovations

Our next generation of economical pumps goes beyond cost.








EN1789 certified for use during transport and E&R scenarios.
Shielded from harsh environments with a validated **IP44** rating.




What makes it new? It's the way you use it.




Upgraded:

-  Bigger screen. Lighter Weight
-  Wider flow rate ranges
-  Higher sensitivity and accuracy
-  Larger storage capacity
-  Longer battery endurance



The New Multifunction Comprehensive Solution



Always Ready to Go.






Functionally Extensive. Clinically Inclusive.

10 different working modes with wider setting ranges:

Rate/Time/Weight/Drip^[1]/Dose Time/Intermittent/Ramp/Micro/First Dose/Sequential

- Safer for pediatric/neonatal use
- Smooth & flexible workflow - up to **9** phases of ramp up/down and up to **10** sequential setups available
- Rates of up to **2000**ml/h
- Infusion Time up to **~100**hrs^[2]




Recommended use for some vasodilators, (eg. Nitroglycerin and triclopidine), to prevent abrupt discontinuation causes rebound vasoconstriction and angina attacks strains.

* [1] Drip mode only applicable for ME660.


[2] Time set up 00:00:01~99:59:59 hh:mm:ss

M260

Enhanced Sensitivity and Accuracy with Strong Functionality.



Robust emergency support with excellent portability.



- Super lightweight: only **1.6kg**^[4] with handle
- Longer battery endurance: up to **12hrs**^[5]
- Adjustment for multi-environment: **10** toggleable levels for screen brightness and alarm volume
- **QuickStart** protocol to boost aid-giving

^[3] When infusion rate ≥ 1.00ml/h

^[4] Including built-in battery

^[5] With optional 5100mAh battery, under rates of 5ml/h


SafeGuard: Risk Mitigation Solutions.

- Drug Library**
Up to 3,000 drug names
- Infusion Logs**
No less than 2,000 records
- Cumulant Management**
4 different auto cumulant modes^[6]

Patented DPS system

Dynamic Pressure Monitoring and Automatic Pressure Release system.

Up to 15 selectable pressure release triggers
Real-time dynamic pressure detection



^[6] Auto cumulant modes include: Recent cumulant/Timed cumulant/Interval cumulant/24h cumulant.

